


## Procedure - Au Exfoliation of TMD

This technique allows us to disassemble 2D materials layer by layer to obtain a monolayer of macroscopic dimension.



The procedure can be split into the following stages:

**Preparatory Stage: Ensure the materials and substrates needed for mechanical exfoliation are ready.**

- ☐ Deposit 150 nm gold<sup>1</sup> layer on an ultra-flat silicon wafer (NOVA Electronic Materials LLC, p-type doped) via e-beam evaporation (Angstrom Engineering EvoVac Multi-Process thin film deposition system) **CONDUCT IN CLEANROOM**
  - Silicon should be roughly 4 inches; O<sub>2</sub> plasma cleaned for 20 minutes
  - The gold is supposed to be an atomically flat gold tape (ultra-flat gold film on a polymer substrate) which is achieved using a template-stripping technique at  $10^{-7}$  torr and 0.5 angstroms per second depositing (this guideline does not need to be as stringently followed since it can take a significant amount of time)
  - **How long can we keep these gold/silicon pieces before they become unusable?**
- ☐ Spin-coat the Au film with a layer of polyvinylpyrrolidone (PVP) (Sigma Aldrich, mw 40000, 10% wt in ethanol/acetonitrile wt 1/1) <https://www.sigmaaldrich.com/US/en/product/sial/pvp40>
  - 3000 rpm, acceleration 1000 rpm/s, 2min
  - Put the silicon with Au and PVP on a hotplate at 150 °C for 5 minutes
  - The point of this PVP layer is to prevent TRT residue from getting on everything
  - 4.5 g ethanol; 4.5 g acetonitrile; 1 g PVP to make the solution

<sup>1</sup> The ultraflat gold tape allows intimate and uniform vdW contact between the gold and the 2D vdW crystal surface, exfoliating a complete monolayer that can be transferred onto the desired substrate. Roughness of gold is not an issue. **What other metals aside from Au have been used for this method?**

### Exfoliation Stage: Manually exfoliate the monolayer.

- ☐ Pick up PVP and gold using thermal release tape (TRT) (semiconductor corp. Release temperature 90 °C); peel everything off of the silicon
  - Gentle pressing to cover the entire area
- ☐ Press the gold side of this layered TRT/PVP/Au onto the surface of a bulk vdW crystal<sup>2</sup>
  - Use tape to pick up a single crystal and fix it onto IPA cleaned glass slide, cut out the rest of the tape so no residue is picked up
  - Apply gentle pressure
  - We want a reasonably sized bulk crystal — too small of a bulk crystal makes the method harder to implement — use a relatively large crystal cleaved with tape (always use an inner surface of the crystal)
  - At this point, you want to complete the procedure within 12 hours since your monolayer has been exposed
- ☐ Peel off the TRT/PVP/Au from the bulk — now there should be a monolayer on the Au. Transfer this to a silicon substrate.
  - Peel slowly as we don't want cracks or folds on the sample
  - **NO** O<sub>2</sub> plasma (UV) cleaning of this new silicon substrate!!! We want the monolayer to stick

### Cleaning Stage: Remove the TRT/PVP/Au layers to expose the monolayer.

- ☐ Place substrate with monolayer onto hotplate at 130 degrees Celsius
- ☐ Place substrate face up in deionized (DI) water to dissolve PVP for 2 hours
- ☐ The sample on the substrate covered by Au is rinsed with acetone (to dry it use gentle N<sub>2</sub> or let the fume hood dry it) and cleaned by O<sub>2</sub> plasma for a little less than 3 minutes (not any longer) to remove any remaining polymer residues.
- ☐ Place substrate face up in etchant solution( I<sub>2</sub>/I<sup>-</sup>) for 10 minutes to dissolve Au
  - KI/I<sub>2</sub> gold etchant solution (2.5g I<sub>2</sub> and 10g KI in 100ml DI water. Iodine, 99.99%, Alfa Aesar; potassium iodide, 99.9%, Alfa Aesar) as mentioned in [1]
  - We do not know if this solution will impact the optical properties of the sample so it is best to not keep the substrate in the etchant longer than 10 minutes
- ☐ Now we should have a monolayer! Rinse it with DI water and isopropanol, and dry it with N<sub>2</sub>.
  - Place sample face up in DI water for 1 hour
  - Quickly rinse in IPA and dry with N<sub>2</sub>

---

<sup>2</sup> Pick up hBN with a different metal. If attempting to use graphene, deposit gold directly onto bulk graphene and pick up the graphene in that way.

## Inventory of Supplies

### Gold Deposition:

- gold crucible
- UV-cleaned silicon wafer (whole)

### Mechanical Exfoliation:

- PVP solution — 4.5 g ethanol (in fume hood; waste container TBD) 4.5 g acetonitrile (need waste container) 1 g pvp powder (chemical bench cabinet 1)
- as of now, there is an existing pvp solution in the fume hood
- thermal release tape (exfoliation drawer)

### Cleaning:

- deionized water (cabinet under fume hood)
- acetone, IPA, N2 gun
- etchant solution — iodine (TBD- need waste container) , potassium iodide (cabinet under fume hood)

A similar procedure for gold-tape exfoliation of TMD from [2]:

- “We first deposit a 150 nm gold film on Si wafer by e-beam evaporation (0.05 nm/s), then spin-coat the gold film by polyvinylpyrrolidone (PVP) solution (Sigma Aldrich, mw 40000, 10% wt in ethanol/acetonitrile wt 1/1) at 1500 rpm for 2 min, with the acceleration of 500 rpm/s, and anneal at 150 °C for 2 min. A single-sided heat release tape is cut into small pieces (~1 × 1 cm<sup>2</sup> square) and stuck onto the PVP/ gold surface to peel off the gold film from the Si wafer to form a gold tape, which is then gently pressed onto a TMD single crystal (purchased from HQ Graphene) to exfoliate a monolayer. The TMD monolayer on the gold tape is then transferred onto a desired substrate. We then heat the substrate with everything on top using a hot plate at 135 °C for 3 min to remove the heat release tape, followed by water soaking for 4 h to remove the PVP layer on gold. Finally, the gold layer is removed by gold etchant (2.5 g I<sub>2</sub> and 10 g KI in 100 mL deionized water), and the TMD monolayer on substrate is washed by water and isopropanol, then dried by a nitrogen gun.”

## Sources

[1] Liu F, Wu W, Bai Y, Chae SH, Li Q, Wang J, Hone J, Zhu XY. Disassembling 2D van der Waals crystals into macroscopic monolayers and reassembling into artificial lattices. *Science*. 2020 Feb 21;367(6480):903-906. doi: 10.1126/science.aba1416. PMID: 32079769. <https://par.nsf.gov/servlets/purl/10155494>

[2] Li, Q., Alfrey, A., Hu, J. *et al.* Macroscopic transition metal dichalcogenides monolayers with uniformly high optical quality. *Nat Commun* 14, 1837 (2023). <https://doi.org/10.1038/s41467-023-37500-1>